

Anti-Prame Polyclonal Antibody For Wb And Elisa - P78395

Item Number: CM0012942

Introduction

High-quality rabbit polyclonal antibody against human PRAME (melanoma antigen). Validated for WB and ELISA; cross-reacts with mouse and rat. Ideal for cancer, apoptosis, and ubiquitination studies.

[Learn More](#)

Parameter	Value
Product Name	PRAME Rabbit pAb
Remarks/Alias	MAPE; OIP4; CT130; OIP-4; PRAME
Species	Human
Gene ID (Human)	23532
Gene ID	23532
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human PRAME (NP_006106.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P78395
Protein Weight	58kDa
Shipping	Ice bag