

Anti-Prb4 Polyclonal Antibody For Wb, Elisa - P10163

Item Number: CM0023955

Introduction

Rabbit polyclonal antibody against human basic salivary proline-rich protein 4 (PRB4, parotid o protein). Suitable for WB and ELISA with cross-reactivity to human and mouse. Ideal for salivary proteomics and oral biology research.

[Learn More](#)

Parameter	Value
Product Name	PRB4 Rabbit pAb
Remarks/Alias	Po; PRB4
Species	Human
Gene ID (Human)	5545
Gene ID	5545
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 178-247 of human PRB4 (NP_002714.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P10163
Protein Weight	31kDa
Shipping	Ice bag