

Anti-Prc1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O43663

Item Number: CM0028441



Introduction

Rabbit polyclonal anti-PRC1 antibody validated for WB, IHC-P, IF/ICC, and ELISA. Recognizes PRC1 in human, mouse, and rat. Suitable for studying cytokinesis, cell division, and bladder cancer research. MW ~72 kDa.

[Learn More](#)

Parameter	Value
Product Name	PRC1 Rabbit pAb
Remarks/Alias	ASE1; PRC1
Species	Human
GeneID (Human)	9055
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 351-620 of human PRC1 (NP_003972.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43663
Protein Weight	72kDa
Shipping	Ice bag