

Anti-Prcc Polyclonal Antibody For Wb And Elisa - Q92733

Item Number: CM0021119

Introduction

Anti-PRCC polyclonal antibody validated for WB and ELISA applications. Reacts with human PRCC protein (UniProt Q92733). Ideal for studying cell cycle regulation and papillary renal cell carcinoma.

[Learn More](#)

Parameter	Value
Product Name	Papillary renal cell carcinoma Rabbit pAb
Remarks/Alias	TPRC; RCCP1; Papillary renal cell carcinoma
Species	Human
GeneID (Human)	5546
GeneID	5546
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 240-491 of human Papillary renal cell carcinoma (NP_005964.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q92733
Protein Weight	52kDa
Shipping	Ice bag