


Anti-Prd Rabbit Monoclonal Antibody For Wb, Elisa - P12955

Item Number: CM0002315



Introduction

High-quality rabbit monoclonal antibody against human PRD (prolidase, PEPD), suitable for WB and ELISA applications, with cross-reactivity to mouse. Ideal for collagen metabolism and prolidase research.

[Learn More](#)

Parameter	Value
Product Name	PRD Rabbit mAb
Remarks/Alias	PROLIDASE; PRD
Species	Human
GeneID (Human)	5184
GeneID	5184
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 400-493 of human PRD (NP_000276.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P12955
Protein Weight	55kDa
Shipping	Ice bag