

Anti-Prdm1 Polyclonal Antibody For Wb, If/Icc, Elisa - 075626

Item Number: CM0023248



Introduction

Rabbit polyclonal antibody targeting PRDM1 (BLIMP-1), validated for western blot, immunofluorescence/ICC, and ELISA. Detects human and mouse PRDM1; ideal for studying transcriptional regulation in immune cell differentiation and tissue-resident lymphocyte biology.

[Learn More](#)

Parameter	Value
Product Name	PRDM1 Rabbit pAb
Remarks/Alias	BLIMP1; PRDI-BF1; PRDM1
Species	Human
Gene ID (Human)	639
Gene ID	639
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 70-210 of human PRDM1 (NP_001189.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O75626
Protein Weight	92kDa
Shipping	Ice bag