

Anti-Prdm11 Polyclonal Antibody For Wb, Elisa - Q9Nqv5

Item Number: CM0029495



Introduction

High-quality rabbit polyclonal antibody targeting human PRDM11 (PFM8), validated for Western blotting and ELISA. Shows cross-reactivity with mouse. Ideal for studying transcription regulation. UniProt: Q9NQV5.

[Learn More](#)

Parameter	Value
Product Name	PRDM11 Rabbit pAb
Remarks/Alias	PFM8; PRDM11
Species	Human
Gene ID	56981
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9NQV5
Protein Weight	58kDa
Shipping	Ice bag