

Anti-Prdm12 Polyclonal Antibody For Wb, Elisa - Q9H4Q4

Item Number: CM0019949



Introduction

PRDM12 Rabbit pAb is a polyclonal antibody for WB and ELISA, detecting human and mouse PRDM12. Essential for nociceptive neuron development research and pain pathway studies.

[Learn More](#)

Parameter	Details
Product Name	PRDM12 Rabbit pAb
Remarks/Alias	PFM9; HSN8; PRDM12
Species	Human
GeneID (Human)	59335
GeneID	59335
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-330 of human PRDM12 (NP_067632.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H4Q4
Protein Weight	40kDa
Shipping	Ice bag