

Anti-Prdm14 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Gzv8

Item Number: CM0014734

Introduction

Rabbit polyclonal antibody targeting human PRDM14, suitable for WB, IHC-P, and ELISA. Detects human, mouse, and rat PRDM14. Ideal for stem cell pluripotency and germ cell development research.

[Learn More](#)

Parameter	Value
Product Name	PRDM14 Rabbit pAb
Remarks/Alias	PFM11; PRDM14
Species	Human
GeneID (Human)	63978
GeneID	63978
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human PRDM14 (NP_078780.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZV8
Protein Weight	64kDa
Shipping	Ice bag