

Anti-Prdm14 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa, Chip - Q9Gzv8

Item Number: CM0030260

Introduction

PRDM14 rabbit polyclonal antibody validated for WB, IHC-P, ELISA, CHIP in human, mouse, rat. Essential tool for studying pluripotency, embryonic stem cell identity, and germ cell development. 64kDa. Q9GZV8.

[Learn More](#)

Parameter	Value
Product Name	PRDM14 Rabbit pAb
Remarks/Alias	PFM11; PRDM14
Species	Human
Gene ID (Human)	63978
Gene ID	63978
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human PRDM14 (NP_078780.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA, CHIP
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZV8
Protein Weight	64kDa
Shipping	Ice bag