

Anti-Prdm6 Polyclonal Antibody For Wb, Elisa - Q9Nqx0

Item Number: CM0029688

Introduction

Anti-PRDM6 rabbit polyclonal antibody validated for WB and ELISA, cross-reactive with human, mouse, and rat. Ideal for studying PRDM6 function in smooth muscle differentiation and chromatin regulation.

[Learn More](#)

Parameter	Value
Product Name	PRDM6 Rabbit pAb
Remarks/Alias	PDA3; KMT8C; PRISM; PRDM6
Species	Human
GeneID (Human)	93166
GeneID	93166
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 376-595 of human PRDM6 (NP_001129711.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NQX0
Protein Weight	64kDa
Shipping	Ice bag