

Anti-Prdm7 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nqw5

Item Number: CM0027460

Introduction

Rabbit polyclonal antibody targeting human PRDM7 (PFM4/ZNF910), raised against recombinant protein (aa 210-310). Validated for WB, IF/ICC, ELISA. Epigenetic research tool for studying histone methyltransferase activity and chromatin remodeling.

[Learn More](#)

Parameter	Value
Product Name	PRDM7 Rabbit pAb
Remarks/Alias	PFM4; ZNF910; PRDM7
Species	Human
GeneID (human)	11105
GeneID	11105
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 210-310 of human PRDM7 (NP_001091643.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9NQW5
Protein Weight	56kDa
Shipping	Ice bag