


Anti-Prdm9 Polyclonal Antibody For If-P, Elisa - Q96Eq9

Item Number: CM0024664



Introduction

Rabbit polyclonal antibody targeting mouse Prdm9, validated for IF-P and ELISA. Cross-reacts with mouse and rat; detects PRDM9, a histone methyltransferase critical for meiotic recombination.

[Learn More](#)

Parameter	Specification
Product Name	Prdm9 Rabbit pAb
Remarks/Alias	Rcr1; Dsbc1; repro7; Meisetz; PRDM9-B; G1-419-29; Prdm9
Species	Mouse
GeneID	213389
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 87-370 of mouse Prdm9 (NP_659058.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q96EQ9
Protein Weight	97kDa
Shipping	Ice bag