


Anti-Preb Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q9Hcu5

Item Number: CM0006317



Introduction

The PREB rabbit monoclonal antibody detects human, mouse, and rat PREB/SEC12 protein. Validated for Western blot, IF/ICC, and ELISA, this antibody targets a 45 kDa GEF critical for ER-Golgi transport and transcriptional regulation. Order from CamelBio for reliable IVD sourcing.

[Learn More](#)

Parameter	Value
Product Name	PREB Rabbit mAb
Remarks/Alias	SEC12; PREB
Species	Human
GeneID (Human)	10113
GeneID	10113
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 89-201 of human PREB (NP_037520.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HCU5
Protein Weight	45kDa
Shipping	Ice bag