

Anti-Prelid1 Polyclonal Antibody For Ihc-P, If/Icc, Elisa - Q9Y255

Item Number: CM0020432



Introduction

Rabbit polyclonal antibody against human PRELID1 (25kDa) for IHC-P, IF/ICC, and ELISA. Involved in mitochondrial cardiolipin accumulation and apoptosis regulation. Validated in human samples.

[Learn More](#)

Parameter	Value
Product Name	PRELID1 Rabbit pAb
Remarks/Alias	PX19; PRELI; SBB112; CGI-106
Species	Human
GeneID (Human)	27166
GeneID	27166
Immunogen	Recombinant protein corresponding to a sequence within amino acids 174-219 of human PRELID1 (NP_037369.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9Y255
Protein Weight	25kDa
Shipping	Ice bag