

Anti-Prelid3A Polyclonal Antibody For Ihc-P, Elisa - Q96N28

Item Number: CM0021475



Introduction

Rabbit polyclonal antibody targeting human PRELID3A (SLMO1), a mitochondrial phosphatidic acid transporter. Validated for IHC-P and ELISA applications to support cardiolipin synthesis research.

[Learn More](#)

Parameter	Value
Product Name	PRELID3A Rabbit pAb
Remarks/Alias	SLMO1; C18orf43; HFL-EDDG1
Species	Human
Gene ID (Human)	10650
Gene ID	10650
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 70-172 of human PRELID3A (NP_006544.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q96N28
Protein Weight	19kDa
Shipping	Ice bag