


Anti-Prelid3B Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y3B1

Item Number: CM0021775



Introduction

Rabbit polyclonal antibody against human PRELID3B (SLMO2). Suitable for Western blot and ELISA. Raised against a synthetic peptide from amino acids 80-180. Useful for studying this PRELI domain-containing phosphoprotein.

[Learn More](#)

Parameter	Value
Product Name	PRELID3B Rabbit pAb
Remarks/Alias	SLMO2, C20orf45, dj543j19.5, PRELID3B
Species	Human
Gene ID (Human)	51012
Gene ID	51012
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 80-180 of human PRELID3B (NP_057129.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y3B1
Protein Weight	21kDa
Shipping	Ice bag