

Anti-Prelp Polyclonal Antibody For Wb And Elisa - P51888

Item Number: CM0022802

Introduction

High-quality PRELP rabbit polyclonal antibody validated for WB and ELISA. Reacts with human, mouse, and rat PRELP. Suitable for extracellular matrix and basement membrane studies.

[Learn More](#)

Parameter	Value
Product Name	PRELP Rabbit pAb
Remarks/Alias	MST161; SLRR2A; MSTP161; PRELP
Species	Human
Gene ID (Human)	5549
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-360 of human PRELP (NP_002716.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P51888
Protein Weight	44kDa
Shipping	Ice bag