


## Anti-Presenilin 1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P49768

Item Number: CM0024207



### Introduction

Rabbit polyclonal antibody against human Presenilin 1, validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for Alzheimer's disease research and gamma-secretase studies.

[Learn More](#)

| Parameter        | Value                                                                                                                     |
|------------------|---------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Presenilin 1 Rabbit pAb                                                                                                   |
| Remarks/Alias    | AD3; FAD; PS1; PS-1; S182; PSNL1; ACNINV3; Presenilin 1                                                                   |
| Species          | Human                                                                                                                     |
| Gene ID (Human)  | 5663                                                                                                                      |
| Gene ID          | 5663                                                                                                                      |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 290-380 of human Presenilin 1 (NP_000012.1) |
| Source           | Rabbit                                                                                                                    |
| Category         | Polyclonal Antibodies                                                                                                     |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                                  |
| Cross Reactivity | Human, Mouse, Rat                                                                                                         |
| SWISS            | P49768                                                                                                                    |
| Protein Weight   | 53kDa                                                                                                                     |
| Shipping         | Ice bag                                                                                                                   |