


Anti-Prex1 Polyclonal Antibody For Wb, Elisa - Q8Tcu6

Item Number: CM0019188



Introduction

PREX1 Rabbit pAb detects human PREX1, a RAC GEF activated by PIP3 and G-protein beta-gamma subunits. Suitable for WB and ELISA; ideal for signal transduction and immunology research.

[Learn More](#)

Parameter	Value
Product Name	PREX1 Rabbit pAb
Remarks/Alias	P-REX1; PREX1
Species	Human
Gene ID (Human)	57580
Immunogen	Recombinant protein containing a sequence corresponding to amino acids 1000-1300 of human PREX1 (NP_065871.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8TCU6
Protein Weight	186kDa
Shipping	Ice bag