

Anti-Prg2 Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P13727

Item Number: CM0028998



Introduction

Anti-PRG2 rabbit polyclonal antibody validated for WB, IF/ICC, IHC-P, and ELISA. Detects human, mouse, and rat PRG2. Ideal for studying immune defense mechanisms, hypersensitivity reactions, and antiparasitic responses.

[Learn More](#)

Product Name	Remarks/Alias	Species	GeneID (Human)	GeneID	Immunogen	Source	Category	Application	Cross Reactivity	SWISS	Protein Weight	Shipping
PRG2 Rabbit pAb	MBP; BMPG; MBP1; proMBP; PRG2	Human	5553	5553	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human PRG2 (NP_002719.3)	Rabbit	Polyclonal Antibodies	WB, IF/ICC, IHC-P, ELISA	Human, Mouse, Rat	P13727	25 kDa	Ice bag