

## Anti-Prg4 Rabbit Polyclonal Antibody For Wb And Elisa - Q92954

Item Number: CM0013563

### Introduction

Rabbit polyclonal antibody targeting human PRG4 (lubricin/MSF). Suitable for Western blot and ELISA, cross-reactive with rat. Predicted MW 151 kDa; useful for studying joint lubrication and synovial biology.

[Learn More](#)

| Parameter        | Value                                                                                                                  |
|------------------|------------------------------------------------------------------------------------------------------------------------|
| Product Name     | PRG4 Rabbit pAb                                                                                                        |
| Remarks/Alias    | MSF; SZP; CACP; HAPO; JCAP; PRG4                                                                                       |
| Species          | Human                                                                                                                  |
| Gene ID          | 10216 (Human)                                                                                                          |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1264-1363 of human PRG4 (NP_001121180.2) |
| Source           | Rabbit                                                                                                                 |
| Category         | Polyclonal Antibodies                                                                                                  |
| Application      | WB, ELISA                                                                                                              |
| Cross Reactivity | Human, Rat                                                                                                             |
| SWISS            | Q92954                                                                                                                 |
| Protein Weight   | 151kDa                                                                                                                 |
| Shipping         | Ice bag                                                                                                                |