

Anti-Prim2 Rabbit Polyclonal Antibody For Wb, Elisa - P49643

Item Number: CM0015133



Introduction

Rabbit polyclonal antibody targeting human PRIM2 (DNA primase large subunit). Validated for Western blot and ELISA. Cross-reacts with human and mouse. Ideal for DNA replication and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	PRIM2 Rabbit pAb
Remarks/Alias	p58; PRIM2A; PRIM2
Species	Human
Gene ID (Human)	5558
Gene ID	5558
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human PRIM2 (NP_000938.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P49643
Protein Weight	59kDa
Shipping	Ice bag