


Anti-Prion Protein (Cd230/Prpc) Rabbit Polyclonal Antibody For Wb, Elisa - P04156

Item Number: CM0028380



Introduction

Rabbit polyclonal antibody targeting human Prion Protein (PrP, CD230), validated for WB and ELISA, with cross-reactivity to mouse. Suitable for prion disease research and diagnostic development.

[Learn More](#)

Parameter	Specification
Product Name	Prion Protein Rabbit pAb
Remarks/Alias	CJD; GSS; PrP; ASCR; KURU; PRIP; PrPc; CD230; AltPrP; p27-30; PrP27-30; PrP33-35C; Prion Protein
Species	Human
GeneID	5621
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-230 of human Prion Protein (NP_001073592.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P04156
Protein Weight	28kDa
Shipping	Ice bag