

Anti-Prion Protein Monoclonal Antibody For Wb, Elisa - P04156

Item Number: CM0008615



Introduction

Rabbit monoclonal antibody targeting human Prion Protein (P04156). Validated for Western blot (WB) and ELISA, cross-reactive with mouse and rat. Suitable for prion disease and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	Prion Protein Rabbit mAb
Remarks/Alias	CJD; GSS; PrP; ASCR; KURU; PRIP; PrPc; CD230; AltPrP; p27-30; PrP27-30; PrP33-35C; Prion Protein
Species	Human
Gene ID (Human)	5621
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 154-253 of human Prion Protein (P04156).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P04156
Protein Weight	28kDa
Shipping	Ice bag