

## Anti-Prion Protein Polyclonal Antibody For Wb, Elisa - P04156

Item Number: CM0024107



### Introduction

Rabbit polyclonal antibody targeting human Prion Protein (PRNP). Validated for Western blot and ELISA, cross-reacts with mouse. Useful for research on prion diseases, neuronal development, and iron homeostasis.

[Learn More](#)

| Parameter        | Value                                                                                                                        |
|------------------|------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Prion Protein Rabbit pAb                                                                                                     |
| Remarks/Alias    | CJD; GSS; PrP; ASCR; KURU; PRIP; PrPc; CD230; AltPrP; p27-30; PrP27-30; PrP33-35C; Prion Protein                             |
| Species          | Human                                                                                                                        |
| GeneID (Human)   | 5621                                                                                                                         |
| GeneID           | 5621                                                                                                                         |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 23-230 of human Prion Protein (NP_001073592.1) |
| Source           | Rabbit                                                                                                                       |
| Category         | Polyclonal Antibodies                                                                                                        |
| Application      | WB, ELISA                                                                                                                    |
| Cross Reactivity | Mouse                                                                                                                        |
| SWISS            | P04156                                                                                                                       |
| Protein Weight   | 28kDa                                                                                                                        |
| Shipping         | Ice bag                                                                                                                      |