


## Anti-Prion Protein Rabbit Polyclonal Antibody For Wb, Elisa - P04156

Item Number: CM0023536



### Introduction

KO-validated rabbit polyclonal antibody against human Prion Protein (PRNP, P04156). Suitable for WB and ELISA, cross-reacts with human, mouse, and rat. For neurodegeneration and prion biology studies.

[Learn More](#)

| Parameter        | Value                                                                                                                         |
|------------------|-------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | [KO Validated] Prion Protein Rabbit pAb                                                                                       |
| Remarks/Alias    | CJD; GSS; PrP; ASCR; KURU; PRIP; PrPc; CD230; AltPrP; p27-30; PrP27-30; PrP33-35C; in                                         |
| Species          | Human                                                                                                                         |
| GeneID (human)   | 5621                                                                                                                          |
| GeneID           | 5621                                                                                                                          |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 23-230 of human Prion Protein (NP_001073592.1). |
| Source           | Rabbit                                                                                                                        |
| Category         | Polyclonal Antibodies                                                                                                         |
| Application      | WB, ELISA                                                                                                                     |
| Cross Reactivity | Human, Mouse, Rat                                                                                                             |
| SWISS            | P04156                                                                                                                        |
| Protein Weight   | 28kDa                                                                                                                         |
| Shipping         | Ice bag                                                                                                                       |