


Anti-Prk2/Pkn2 Monoclonal Antibody For Wb And Elisa - Q16513

Item Number: CM0002768



Introduction

Rabbit monoclonal antibody against human PRK2/PKN2 (PKN gamma), suitable for WB and ELISA. Cross-reacts with mouse and rat. Ideal for research in cell cycle, cytoskeleton dynamics, and cell migration.

[Learn More](#)

Specification	Details
Product Name	PRK2/PKN2 Rabbit mAb
Remarks/Alias	PAK2; PRK2; STK7; Pak-2; PRKCL2; PRO2042; PRK2/PKN2
Species	Human
GeneID (Human)	5586
GeneID	5586
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16513
Protein Weight	112kDa
Shipping	Ice bag