


Anti-Prk2/Pkn2 Rabbit Polyclonal Antibody For Wb And Elisa - Q16513

Item Number: CM0022209



Introduction

Rabbit polyclonal antibody targeting human PRK2/PKN2, validated for WB and ELISA. Cross-reacts with mouse, rat. Suitable for investigating Rho/Rac signaling, cell cycle progression, and actin cytoskeleton dynamics.

[Learn More](#)

Parameter	Value
Product Name	PRK2/PKN2 Rabbit pAb
Remarks/Alias	PAK2; PRK2; STK7; Pak-2; PRKCL2; PRO2042; PRK2/PKN2
Species	Human
GeneID (Human)	5586
GeneID	5586
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q16513
Protein Weight	112kDa
Shipping	Ice bag