

Anti-Prkacg Rabbit Polyclonal Antibody For Wb & Elisa - P22612

Item Number: CM0024888

Introduction

High-quality rabbit polyclonal antibody targeting human PRKACG for Western blotting and ELISA. Cross-reacts with mouse and rat. Ideal for cAMP-dependent kinase signaling and testis research.

[Learn More](#)

Parameter	Value
Product Name	PRKACG Rabbit pAb
Remarks/Alias	KAPG; PKACg; BDPLT19; PRKACG
Species	Human
Gene ID (Human)	5568
Gene ID	5568
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-351 of human PRKACG (NP_002723.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22612
Protein Weight	40kDa
Shipping	Ice bag