


Anti-Prkag1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P54619

Item Number: CM0028303



Introduction

Rabbit polyclonal antibody targeting human PRKAG1 (AMPK gamma-1 subunit). Validated for WB, IHC-P, IF/ICC, ELISA with human, mouse, rat reactivity. Suitable for energy metabolism and cancer research.

[Learn More](#)

Parameter	Value
Product Name	PRKAG1 Rabbit pAb
Remarks/Alias	AMPKG; PRKAG1
Species	Human
Gene ID (Human)	5571
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-331 of human PRKAG1 (NP_002724.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54619
Protein Weight	38kDa
Shipping	Ice bag