

Anti-Prkag3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Ugi9

Item Number: CM0017987

Introduction

Rabbit polyclonal antibody targeting human PRKAG3 (AMPK gamma-3 subunit) for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studying cellular energy metabolism and muscle energy regulation.

[Learn More](#)

Parameter	Value
Product Name	PRKAG3 Rabbit pAb
Remarks/Alias	AMPKG3; SMGMQTL; PRKAG3
Species	Human
GeneID (Human)	53632
GeneID	53632
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human PRKAG3 (NP_059127.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UGI9
Protein Weight	54kDa
Shipping	Ice bag