

Anti-Prkar2B Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P31323

Item Number: CM0002559

Introduction

Rabbit monoclonal antibody against PRKAR2B, the regulatory subunit of cAMP-dependent protein kinase type II-beta. Validated for WB, IHC-P, ELISA with human, mouse, rat cross-reactivity. Ideal for cAMP signaling studies.

[Learn More](#)

Parameter	Value
Product Name	PRKAR2B Rabbit mAb
Remarks/Alias	PRKAR2; RII-BETA; PRKAR2B
Species	Human
GeneID (Human)	5577
GeneID	5577
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 319-418 of human PRKAR2B (P31323).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P31323
Protein Weight	46kDa
Shipping	Ice bag