

Anti-Prkcdbp (Cavin-3) Rabbit Polyclonal Antibody For Wb, Elisa - Q969G5

Item Number: CM0020684

Introduction

High-quality rabbit polyclonal antibody targeting human PRKCDBP (Cavin-3/SRBC), validated for Western blot and ELISA. Cross-reacts with rat. Suitable for caveolae dynamics and circadian rhythm research.

[Learn More](#)

Parameter	Value
Product Name	PRKCDBP Rabbit pAb
Remarks/Alias	SRBC; HSRBC; PRKCDBP; cavin-3
Species	Human
GeneID (Human)	112464
Immunogen	Recombinant Protein: Recombinant fusion protein containing a sequence corresponding to amino acids 75-160 of human PRKCDBP (NP_659477.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q969G5
Protein Weight	28kDa
Shipping	Ice bag