

Anti-Prkcsh Monoclonal Antibody For Wb, Elisa - P14314

Item Number: CM0003818



Introduction

Anti-PRKCSH rabbit monoclonal antibody validated for WB and ELISA, specifically reacting with human PRKCSH (59kDa), a regulatory subunit of glucosidase II in the endoplasmic reticulum. Ideal for research on glycoprotein quality control and ER functions.

[Learn More](#)

Parameter	Value
Product Name	PRKCSH Rabbit mAb
Remarks/Alias	GIIB; PCLD; PLD1; G19P1; PCLD1; PKCSH; AGE-R2; VASAP-60; PRKCSH
Species	Human
GeneID (Human)	5589
GeneID	5589
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 429-528 of human PRKCSH (P14314).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P14314
Protein Weight	59kDa
Shipping	Ice bag