


Anti-Prkd3 Rabbit Polyclonal Antibody For Wb And Elisa - 094806

Item Number: CM0028057



Introduction

Rabbit polyclonal antibody targeting human PRKD3 (serine/threonine-protein kinase D3, EPK2, PRKCN) for WB and ELISA applications. Cross-reacts with mouse. Ideal for diacylglycerol signaling and oxidative stress research.

[Learn More](#)

Parameter	Value
Product Name	PRKD3 Rabbit pAb
Remarks/Alias	EPK2; PKD3; PRKCN; PKC-NU; nPKC-NU; PRKD3
Species	Human
GeneID	23683
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human PRKD3 (NP_005804.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	094806
Protein Weight	100kDa
Shipping	Ice bag