


Anti-Prkg1 Polyclonal Antibody For Wb And Elisa - Q13976

Item Number: CM0028855



Introduction

Rabbit polyclonal antibody against human PRKG1 (cGMP-dependent protein kinase 1), suitable for Western blot and ELISA. Cross-reacts with mouse. Essential for NO/cGMP signaling research in smooth muscle, platelet, and CNS studies.

[Learn More](#)

Parameter	Value
Product Name	PRKG1 Rabbit pAb
Remarks/Alias	PKG; cGK; AAT8; PKG1; cGK1; cGKI; cGK 1; PRKG1B; PRKGR1B; cGKI-BETA; cGKI-alpha; PRKG1
Species	Human
Gene ID (Human)	5592
Gene ID	5592
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human PRKG1 (NP_006249.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13976
Protein Weight	76kDa
Shipping	Ice bag