

Anti-Prkg2 Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q13237

Item Number: CM0020802



Introduction

Rabbit polyclonal antibody against human PRKG2 (cGKII) validated for WB, IF/ICC, IP, ELISA. Cross-reacts with mouse and rat. Suited for studies of intestinal secretion, bone growth, and synaptic plasticity.

[Learn More](#)

Parameter	Value
Product Name	PRKG2 Rabbit pAb
Remarks/Alias	AMD4; PKG2; SMDP; cGK2; cGKII; PRKGR2; PRKG2
Species	Human
Gene ID (Human)	5593
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human PRKG2 (NP_006250.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q13237
Protein Weight	87kDa
Shipping	Ice bag