

Anti-Prkra/Pact Rabbit Polyclonal Antibody For Wb, Elisa - 075569

Item Number: CM0028770



Introduction

Rabbit polyclonal antibody targeting human PRKRA/PACT, validated for Western blot (WB) and ELISA. Detects a ~34kDa band and cross-reacts with mouse. Suitable for studying PKR activation, RNA interference, and apoptosis pathways.

[Learn More](#)

Parameter	Value
Product Name	PRKRA/PACT Rabbit pAb
Remarks/Alias	RAX; PACT; DYT16; HSD14; PRKRA/PACT
Species	Human
GeneID (Human)	8575
GeneID	8575
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-313 of human PRKRA/PACT (NP_003681.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	075569
Protein Weight	34kDa
Shipping	Ice bag