


Anti-Prkrip1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H875

Item Number: CM0021636



Introduction

PRKRIP1 Rabbit polyclonal antibody detects human, mouse, and rat PRKRIP1, validated for WB, IF/ICC, and ELISA. PRKRIP1 is a crucial spliceosomal component required for pre-mRNA splicing, making this antibody ideal for RNA processing and spliceosome research.

[Learn More](#)

Parameter	Value
Product Name	PRKRIP1 Rabbit pAb
Remarks/Alias	C114; KRBOX3; PRKRIP1
Species	Human
GeneID (Human)	79706
GeneID	79706
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 67-134 of human PRKRIP1 (NP_078929.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H875
Protein Weight	21kDa
Shipping	Ice bag