

Anti-Prl Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P01236

Item Number: CM0021436



Introduction

Rabbit polyclonal antibody targeting human prolactin (PRL, ~26kDa). Validated for WB, IHC-P, and ELISA; cross-reacts with mouse and rat. Ideal for lactation and mammary gland research.

[Learn More](#)

Parameter	Value
Product Name	PRL Rabbit pAb
Remarks/Alias	GHA1; PRL
Species	Human
Gene ID (Human)	5617
Gene ID	5617
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-227 of human PRL (NP_001157030.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01236
Protein Weight	26kDa
Shipping	Ice bag