

Anti-Prl/Prolactin Rabbit Monoclonal Antibody For Ihc-P, Elisa - P01236

Item Number: CM0009900



Introduction

Recombinant rabbit monoclonal antibody targeting human Prolactin for IHC-P and ELISA applications. Suitable for lactation and hormone research in human samples.

[Learn More](#)

Parameter	Value
Product Name	PRL/Prolactin Rabbit mAb
Remarks/Alias	GHA1
Species	Human
GeneID (Human)	5617
GeneID	5617
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	P01236
Protein Weight	26kDa
Shipping	Ice bag