

Anti-Prmt1 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q99873

Item Number: CM0005347



Introduction

Rabbit monoclonal anti-PRMT1 antibody validated for WB, IHC-P, IF/ICC, ELISA. Recognizes human, mouse, and rat PRMT1 (42kDa). Useful for studying arginine methylation, epigenetic regulation, and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	PRMT1 Rabbit mAb
Remarks/Alias	ANM1; HCP1; IR1B4; HRMT1L2; PRMT1
Species	Human
Gene ID (Human)	3276
Gene ID	3276
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PRMT1 (Q99873).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99873
Protein Weight	42kDa
Shipping	Ice bag