

Anti-Prmt2 Rabbit Polyclonal Antibody For Wb, Elisa - P55345

Item Number: CM0023644



Introduction

High-quality rabbit polyclonal antibody targeting human PRMT2, suitable for Western blot and ELISA applications. Ideal for research in epigenetics and signal transduction involving PRMT2-mediated methylation.

[Learn More](#)

Parameter	Value
Product Name	PRMT2 Rabbit pAb
Remarks/Alias	HRMT1L1; PRMT2
Species	Human
Gene ID	3275
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 290-433 of human PRMT2 (NP_001526.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P55345
Protein Weight	49kDa
Shipping	Ice bag