


Anti-Prmt5 Monoclonal Antibody For Wb, Ihc-P, Elisa - O14744

Item Number: CM0003919



Introduction

High-quality rabbit monoclonal antibody against human PRMT5, suitable for WB, IHC-P, and ELISA. Detects human, mouse, and rat PRMT5 (predicted 73kDa). Ideal for epigenetics and transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	PRMT5 Rabbit mAb
Remarks/Alias	HSL7; JBP1; SKB1; IBP72; SKB1Hs; HRMT1L5; PRMT5
Species	Human
GeneID (Human)	10419
GeneID	10419
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PRMT5 (O14744)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14744
Protein Weight	73kDa
Shipping	Ice bag