


## Anti-Prmt6 Polyclonal Antibody For Wb, If/Icc, Elisa - Q96La8

Item Number: CM0030492



### Introduction

High-quality rabbit polyclonal antibody targeting human PRMT6 (HRMT1L6), suitable for WB, IF/ICC, and ELISA. Cross-reacts with mouse. Supports investigations into histone methylation, transcriptional repression, and DNA damage response.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | [KO Validated] PRMT6 Rabbit pAb                                                                                  |
| Remarks/Alias    | HRMT1L6; T6                                                                                                      |
| Species          | Human                                                                                                            |
| Gene ID (Human)  | 55170                                                                                                            |
| Gene ID          | 55170                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-375 of human PRMT6 (NP_060607.2) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, IF/ICC, ELISA                                                                                                |
| Cross Reactivity | Human, Mouse                                                                                                     |
| SWISS            | Q96LA8                                                                                                           |
| Protein Weight   | 42kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |