

Anti-Prmt6 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q96La8

Item Number: CM0021237



Introduction

PRMT6 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Detects human PRMT6, a histone arginine methyltransferase involved in transcriptional repression, DNA repair, and RNA splicing. Ideal for epigenetics and cancer research.

[Learn More](#)

Parameter	Value
Product Name	PRMT6 Rabbit pAb
Remarks/Alias	HRMT1L6; PRMT6
Species	Human
GeneID (human)	55170
GeneID	55170
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-375 of human PRMT6 (NP_060607.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q96LA8
Protein Weight	42kDa
Shipping	Ice bag