

Anti-Prmt7 Polyclonal Antibody For Wb, Elisa - Q9Nvm4

Item Number: CM0015213

Introduction

Rabbit polyclonal antibody targeting human PRMT7 (Q9NVM4). Validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Suitable for studying histone methylation, snRNP assembly, and gene imprinting.

[Learn More](#)

Parameter	Value
Product Name	PRMT7 Rabbit pAb
Remarks/Alias	SBIDDS; PRMT7
Species	Human
Gene ID	54496 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 513-692 of human PRMT7 (NP_061896.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NVM4
Protein Weight	78kDa
Shipping	Ice bag