

## Anti-Prmt7 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Nvm4

Item Number: CM0022305

### Introduction

Rabbit polyclonal antibody against human PRMT7 (protein arginine N-methyltransferase 7) for WB, IHC-P, ELISA. Research chromatin regulation, gene imprinting, and snRNP assembly. Validated on human samples. 78kDa.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | PRMT7 Rabbit pAb                                                                                                   |
| Remarks/Alias    | SBIDDS; PRMT7                                                                                                      |
| Species          | Human                                                                                                              |
| Gene ID (Human)  | 54496                                                                                                              |
| Gene ID          | 54496                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 513-692 of human PRMT7 (NP_061896.1) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IHC-P, ELISA                                                                                                   |
| Cross Reactivity | Human                                                                                                              |
| SWISS            | Q9NVM4                                                                                                             |
| Protein Weight   | 78kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |