

## Anti-Prmt7 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Nvm4

Item Number: CM0003011

### Introduction

Rabbit monoclonal antibody against human PRMT7 (Q9NVM4). Validated for Western blot and ELISA, cross-reactive with human, mouse, and rat. Ideal for studying arginine methylation and chromatin regulation.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | PRMT7 Rabbit mAb                                                                                                    |
| Remarks/Alias    | SBIDDS; PRMT7                                                                                                       |
| Species          | Human                                                                                                               |
| Gene ID (Human)  | 54496                                                                                                               |
| Gene ID          | 54496                                                                                                               |
| Immunogen        | Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PRMT7 (Q9NVM4). |
| Source           | Rabbit                                                                                                              |
| Category         | Monoclonal Antibodies                                                                                               |
| Application      | WB, ELISA                                                                                                           |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | Q9NVM4                                                                                                              |
| Protein Weight   | 78kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |